

Pain free, Worry free & Hassle free!

The complete Ultra Suction kit

Reduced Acrylic Palate

UltraSuction mounted on an upper denture will help reduce the size of the acrylic palate, improving both phonetics and taste. Clients, no longer embarrassed by loose fitting dentures, may become more confident, eat healthier, be more social and enjoy a better lifestyle.

Regular Patient Recall

UltraSuction helps you establish an ongoing Patient/Practitioner relationship by way of regular recalls for preventative maintenance.

Palateless denture possible



Ultra Suction™
holds your denture



Valve

Processing cap

service key
(gauge)

- Accessories:**
- * A specially designed profiled bar (spacer) used to create a suction chamber along the entire denture at the ridge level.
 - * 2 UltraSuction valves
 - * 2 Processing caps
 - ** Service key
 - ** Spare diaphragms
 - ** Instruction sheet

What is UltraSuction?

Denture wearers can experience an insecure feeling while eating, drinking and talking as their dentures move sideways, backwards and forwards. UltraSuction allows clients a more comfortable experience wearing dentures while enjoying a more sociable lifestyle and building confidence.

UltraSuction is a patented system that hold upper and/or lower dentures in place using a mounted valve and suction chamber. The UltraSuction system comprises denture valves embedded into the denture and connected via two tiny 1mm air passages open to the air chamber. Therefore, there is no contact between the valves and the patient's gum tissue.

Can be built into a NEW denture or into older dentures during the rebase or reline procedure.

How UltraSuction Works.....

UltraSuction works on a fairly simple principle suction. Once the denture is in place and tight, the gum rests against the suction chamber, expelling the air and forming a tight seal. The result - a better fit to the gum tissues and an improved resistance to dislodging forces.



Ultra Suction™
holds your denture

Loose dentures a problem?
UltraSuction holds dentures
firm and secure, worry free!



Improved Speech and Taste More!

UltraSuction mounted on an upper denture helps reduce the size of acrylic palate, thus improving both phonetics and taste.



1-800-767-5303

Visit us on the web @ www.mitchlabs.com

See your
patients smile...
and more often!

